

# IBISXA

**Multifunctional  
manual/foot-propelling  
wheelchair**

*Maximum comfort  
Perfect driving and  
foot-propelling  
30° tilt-in-space  
Great independence*



**Ibis XA. Perfect combination of mobility and comfort for greater independence.** For the semi-active wheelchair user, a smooth and comfortable ride is crucial, and that is exactly what the Ibis XA is offering. This multifunctional manual/foot-propelling wheelchair encourages the individual to make maximum use of his/her own abilities thus enhancing independence. The adjustable wheel position coupled with unparalleled seat width adjustment ensures perfect driving characteristics for every use. Additionally, a low seat height from 39 cm allows easy footpropelling. Our individually adjustable Sedeo seating system offers maximum support and comfort. To ensure the maximum variation in seating posture, the Ibis XA incorporates an easy to operate tilt-in-space adjustment. This can be operated by either the person in the wheelchair or the carer. Fully adjustable, the Ibis XA suits almost everyone and is easy to re-use. Ideal for home and hospital use.



**handicare**

# IBISXA



## Maximum comfort

The fully individually adjustable Sedeo seating system provides excellent physical support and maximum comfort.

## Perfect driving and foot-propelling

The adjustable wheel position coupled with unparalleled seat width adjustment ensures a smooth ride and maximum manoeuvrability. Low seat height (from 39 cm) and plenty of leg space for easy foot-propelling.

## 30° tilt-in-space

The easy to operate tilt-in-space adjustment can be operated by the person in the wheelchair or the carer. With up to 30° tilt-in-space adjustment, optimal variation in seating posture is guaranteed.

## Great independence

The person in the wheelchair can move the wheelchair with great ease – either manually or foot-propelling – thus enhancing freedom and independence.

## DIMENSIONS

Maximum user weight	kg	stone	160	25
Total length	mm	inch	1040	40.9
Total width	mm	inch	660 / 700 /	24.6
Total weight	kg	stone	42	6.3
Seat tilt adjustment	°	°	30	30
Seat angle	°	°	0 - 6	0 - 6
Seat depth	mm	inch	440 - 520	17.3 - 20.5
Seat width	mm	inch	380 - 555	15 - 21.9
Seat height incl. cushion	mm	inch	390 - 540	15.4 - 21.3
Seat height excl. cushion	mm	inch	320 - 470	12.6 - 18.5
Backrest angle	°	°	89 - 120	89 - 120
Back height	mm	inch	520 - 570	20.5 - 22.5
Diameter front wheels	inch	inch	7	7
Diameter rear wheels	inch	inch	20 / 22 / 24	20 / 22 / 24
<b>PERFORMANCE</b>				
Ground clearance	mm	inch	70	2.8
Turning radius (ISO 7176-5)	mm	inch	810	31.9

*The Ibis XA is part of the Ibis X-Series, a range that extends from a push wheelchair, a manual/foot-propelling version to an electric model. All this is based on one, width-adjustable frame with tilt-in-space.*

ISO 7176-19  
EN12182 (1999)  
EN12183 (1999)



handicare

The Netherlands  
T +31 (0)492 593 888  
F +31 (0)492 537 931  
info@handicare.com  
www.handicare.com

UK  
T +44 (0)1932 858 687  
F +44 (0)1932 848 384  
info@handicare.co.uk  
www.handicare.co.uk