

Order form

XP JUNIOR

RETAIL 2012



handicare

T: +44 (0)8450 745 945

F: +44 (0)8450 745 946

mobilityorders@handicare.co.uk

www.handicare.co.uk

Valid from	1-1-2012	Version	1.0
Dealer	Date		
Reference No.			
Delivery address	Quantity		
	PRICE LIST ID	0092	

Standard configuration exclusive VAT £ 4.295,00

- 52 cm / 20.5" seat height
- 28 - 38 cm / 11 - 15" stepless seat width between hip cushions
- 25 - 40 cm / 9.8 - 15.7" stepless seat depth
- Baseframe ice blue
- Rear wheel drive
- Air front & rear wheels
- Shark controller, swing-away bracket
- Optimum Support System, stretch Wales black
- Thoracic supports, stretch Wales black
- Standard swing-away legrests
- Anti-tip brackets
- Approx. 31 Ah batteries
- 5 Ah charger
- Seat belt
- Full range of adjustment options available: seat width - seat depth - seat angle - seat height - centre of gravity - backrest angle.

Colour options can be found in section adaptations or ask your sales representative for other available colours and prices.

The standard configuration below is indicated with 'basic'. If an alternative option is required please mark clearly by ticking the associated box.

THE Y-NUMBERS WHICH ARE MENTIONED BEHIND THE OPTIONS, ARE NOT ARTICLE NUMBERS. ALWAYS USE THE Y-NUMBERS IN COMBINATION WITH THE UNIT-ID NUMBER MENTIONED ABOVE.

BASE FRAME



	Option	Excl.
<input type="checkbox"/> SEAT HEIGHT, INCL. STANDARD CUSHION (8 CM / 3.1"), USER WEIGHT EXCL.		
<input type="checkbox"/> 47,5 cm / 18.7"	Y00573	0
<input type="checkbox"/> 49 cm / 19.3"	Y00186	0
<input type="checkbox"/> 50,5 cm / 19.9"	Y00574	0
<input type="checkbox"/> 52 cm / 20.5"	Y00190	Basic
<input type="checkbox"/> 53,5 cm / 21.1"	Y00193	0
<input type="checkbox"/> 55 cm / 21.7"	Y00196	0
<input type="checkbox"/> 56,5 cm / 22.2"	Y00571	0
<input type="checkbox"/> 58 cm / 22.8"	Y00572	0
<input type="checkbox"/> 1 Rear wheel drive	Y00106	Basic
<input type="checkbox"/> Colour frame ice blue	Y00507	Basic
<input type="checkbox"/> Colour frame Candy Pink	Y00508	0
<input type="checkbox"/> Castors air	Y00097	Basic
<input type="checkbox"/> Castors puncture proof with solid inner tube	Y00098	65
<input type="checkbox"/> Drive wheels air	Y00107	Basic
<input type="checkbox"/> Drive wheels puncture proof with solid inner tube	Y00108	95

ADJUSTMENTS & LIGHTING



	Option	Excl.
<input type="checkbox"/> 1 Electrical backrest adjustment 32°	Y00385	799
<input type="checkbox"/> 2 Electrical tilt adjustment, SH will influence tilt range	Y00398	499
<input type="checkbox"/> 3 Electrical high/low adjustment, min. SH 51 cm / 20.1 (i.c.w. tilt adjustment, see adaptations)	Y00391	1.495
<input type="checkbox"/> 4 Lighting and indicators compact version	Y00395	275



The product picture shown may deviate from the standard configuration.

ISO 7176-19

EN12184 (1999)



CLASS 2

P.O.R. = Price on request, D.O.R. = Delivery time on request

Prices in UK Pounds excluding VAT. This price list will replace all previous price list versions. We reserve the right to change (price) modifications and incorrections.

CONTROLLERS

Option

Excl.




Shark drive only: no lighting or electrical adjustments possible
 Shark lights: lighting, no electrical adjustments possible
 Shark seats & lights: lighting and two electrical adjustments possible
 DX G90A: lighting and five electrical adjustments possible

<input type="checkbox"/>	1 Shark drive only, right	Y00228	Basic
<input type="checkbox"/>	1 Shark drive only, left	Y00227	0
<input type="checkbox"/>	2 Shark lights, right	Y00230	120
<input type="checkbox"/>	2 Shark lights, left	Y00229	120
<input type="checkbox"/>	3 Shark seat & lights, right	Y00232	175
<input type="checkbox"/>	3 Shark seat & lights, left	Y00231	175
<input type="checkbox"/>	4 DX G90A joystick module, right	Y00241	250
<input type="checkbox"/>	4 DX G90A joystick module, left	Y00240	250
<input type="checkbox"/>	5 Attendant control / dual joystick (not i.c.w. shark drive only), right	Y00222	440
<input type="checkbox"/>	5 Attendant control / dual joystick (not i.c.w. shark drive only), left	Y00221	440
<input type="checkbox"/>	6 Swing-away bracket	Y00265	Basic
<input type="checkbox"/>	7 Mushroom knob	Y00267	42

SEAT CUSHIONS

Option

Excl.




<input type="checkbox"/>	Seat cushion width 30 cm / 11.8" for SW 25-33 cm (9.8-13")	Y00489	0
<input type="checkbox"/>	Seat cushion width 35 cm / 13.8" for SW 30-38 cm (11.8-15")	Y00490	Basic
<input type="checkbox"/>	1 Form seat OSS, pre shaped	Y00474	Basic
<input type="checkbox"/>	2 Standard seat	Y00459	0
<input type="checkbox"/>	3 Pan seat, lip height 3 cm / 1.2"	Y00467	60
<input type="checkbox"/>	4 Board seat, 1 cm / 0.4" padded	Y00471	30
<input type="checkbox"/>	5 Webbed seat, extra soft, stretch rubber upholstery	Y00472	150
<input type="checkbox"/>	6 Memory foam, flat, stretch rubber upholstery	Y00461	155
<input type="checkbox"/>	7 Abduction support detachable	Y00442	170

BACKRESTS

Option

Excl.



<input type="checkbox"/>	Backrest cushion width 30 cm / 11.8"	Y00370	0
<input type="checkbox"/>	Backrest cushion width 35 cm / 13.8"	Y00371	Basic
<input type="checkbox"/>	1 Formback OSS, thoracic support included	Y00361	Basic
<input type="checkbox"/>	2 Standard backrest	Y00347	0
<input type="checkbox"/>	3 Webbed backrest, extra soft, stretch fabric upholstery	Y00349	185
<input type="checkbox"/>	4 Backrest flat upholstered	Y00352	25
<input type="checkbox"/>	5 Backrest with adjustable straps, stretch fabric upholstery	Y00357	250
<input type="checkbox"/>	6 Thoracic support on frame (not i.c.w. Formback OSS), right	Y00384	105
<input type="checkbox"/>	6 Thoracic support on frame (not i.c.w. Formback OSS), left	Y00383	105

LEGRESTS

Option

Excl.



<input type="checkbox"/>	COMPLETE SET		
<input type="checkbox"/>	1 Standard swing-away legrests	Y00321	Basic
<input type="checkbox"/>	2 Comfort swing-away legrests with calf plates	Y00334	375

ARMRESTS

Option

Excl.




<input type="checkbox"/>	1 Padded armrest, right	Y00070	Basic
<input type="checkbox"/>	1 Padded armrest, left	Y00057	Basic

SIDE CUSHIONS

Option

Excl.



<input type="checkbox"/>	1 Side cushion padded, stretch Wales, right	Y00438	Basic
<input type="checkbox"/>	1 Side cushion padded, stretch Wales, left	Y00431	Basic
<input type="checkbox"/>	2 Acrylic side cushion, right	Y00437	0
<input type="checkbox"/>	2 Acrylic side cushion, left	Y00430	0
<input type="checkbox"/>	3 Side cushion padded, right	Y00439	5
<input type="checkbox"/>	3 Side cushion padded, left	Y00432	5

HEADRESTS

Option

Excl.



<input type="checkbox"/>	1 Adjustable metal bracket	Y00342	85
<input type="checkbox"/>	1 Width adjustable metal bracket	Y00343	95
<input type="checkbox"/>	2 Headrest cushion small	Y00339	85
<input type="checkbox"/>	3 Headrest cushion flexible	Y00341	90

WORKTOP

Option

Excl.



<input type="checkbox"/>	ENTRY TUBE SIDE NOT I.C.W. REMOTE SIDE		
<input type="checkbox"/>	1 Folding acrylic for joystick recess on the side, entry tube right	Y00415	240
<input type="checkbox"/>	1 Folding acrylic for joystick recess on the side, entry tube left	Y00414	240

