

# Order form

# PUMA 4

## RETAIL 2012



handicare

T: +44 (0)8450 745 945

F: +44 (0)8450 745 946

mobilityorders@handicare.co.uk

www.handicare.co.uk

Valid from	1-1-2012	Version	1.0
Dealer	Date		
Reference No.			
Delivery address	Quantity		
	PRICE LIST ID	0087	

**Standard configuration exclusive VAT £ 4.295,00**

- 49 cm / 19.3" seat height
- 46 cm / 18.1" seat width
- 50 cm / 19.7" seat depth
- Rear wheel drive
- Grey shrouds
- Indoor castors
- Shark controller, fixed bracket
- Standard seat cushion with grey middle section
- Standard backrest cushion with grey middle section
- Standard swing-away legrests with calf strap
- Polyurethane armrests
- Polyurethane side cushions
- Approx. 60 Ah batteries
- 8 Ah charger
- Transport brackets (for hooks & carbiners)
- Seat belt
- Full range of adjustment options available: seat width - seat depth - seat angle - seat height - centre of gravity - backrest angle.



The product picture shown may deviate from the standard configuration.

ISO 7176-19 EN12184 (1999)









CLASS 2

In its standard configuration the chair conforms to the Highways regulations 1988, Class 2 vehicle requirements however it must be noted that when choosing either of the following options, Electrical High/low or Tilt in Space the chair will automatically fall into the Class 3 specification and as such if being used outside the home environment will need the addition of lights indicators and a rear view mirror, this information also applies if you consider adding other accessories that should then make the chair heavier than the 113.4 Kg of a Class 2 Chair.

Colour options can be found in section adaptations or ask your sales representative for other available colours and prices. The standard configuration below is indicated with 'basic'. If an alternative option is required please mark clearly by ticking the associated box.

**THE Y-NUMBERS WHICH ARE MENTIONED BEHIND THE OPTIONS, ARE NOT ARTICLE NUMBERS. ALWAYS USE THE Y-NUMBERS IN COMBINATION WITH THE UNIT-ID NUMBER MENTIONED ABOVE.**





### BASE FRAME

BASE FRAME		Option	Excl.
			
			
<input type="checkbox"/>	SEAT HEIGHT, INCL. STANDARD CUSHION (7 CM / 2.8"), USER WEIGHT EXCL.		
<input type="checkbox"/>	<b>49 cm / 19.3"</b>	<b>Y00186</b>	<b>Basic</b>
<input type="checkbox"/>	51 cm / 20.1"	Y00188	0
<input type="checkbox"/>	53 cm / 20.9"	Y00192	0
<input type="checkbox"/>	55 cm / 21.7"	Y00196	0
<input type="checkbox"/>	1 Front wheel drive	Y00105	0
<input type="checkbox"/>	<b>2 Rear wheel drive</b>	<b>Y00106</b>	<b>Basic</b>
<input type="checkbox"/>	<b>3 Colour shrouds grey</b>	<b>Y00209</b>	<b>Basic</b>
<input type="checkbox"/>	4 Colour shrouds yellow	Y00210	0
<input type="checkbox"/>	Castors type outdoor	Y00095	35
<input type="checkbox"/>	<b>Castors type indoor</b>	<b>Y00096</b>	<b>Basic</b>
<input type="checkbox"/>	<b>Speed 6 km/h / 4 mph</b>	<b>Y00211</b>	<b>Basic</b>
<input type="checkbox"/>	5 Suspended castors for extra stability and comfort (only i.c.w. castors type "outdoor")	Y00219	215
<input type="checkbox"/>	<b>Castors air</b>	<b>Y00097</b>	<b>Basic</b>
<input type="checkbox"/>	Castors puncture proof with solid inner tube (only i.c.w. castors type "outdoor")	Y00098	80
<input type="checkbox"/>	<b>Drive wheels air</b>	<b>Y00107</b>	<b>Basic</b>
<input type="checkbox"/>	Drive wheels puncture proof with solid inner tube	Y00108	100
<input type="checkbox"/>	6 Telescopic Safety Device (only on FWD, obligatory i.c.w. 10 km/h / 6.2 mph)	Y00220	140











P.O.R. = Price on request, D.O.R. = Delivery time on request

Prices in UK Pounds excluding VAT. This price list will replace all previous price list versions. We reserve the right to change (price) modifications and incorrections.










**ADJUSTMENTS & LIGHTING**

		Option	Excl.
	<input type="checkbox"/> 1	Electrical backrest adjustment 32°	Y00385 699
	<input type="checkbox"/> 2	Electrical tilt adjustment 40° SH + 7 cm	Y00398 650
	<input type="checkbox"/> 3	Electrical high/low adjustment 31 cm / 12.2" SH + 7 cm	Y00391 950
	<input type="checkbox"/> 4	Lighting and indicators compact version	Y00395 275
			

**CONTROLLERS**

		Option	Excl.
		Shark drive only: no lighting or electrical adjustments possible	
		Shark lights: lighting, no electrical adjustments possible	
		Shark seats & lights: lighting and two electrical adjustments possible	
		DX G90A: lighting and five electrical adjustments possible	
		DX 12-key flat: lighting and five electrical adjustments possible	
	<input type="checkbox"/> 1	<b>Shark drive only, right</b>	<b>Y00228 Basic</b>
	<input type="checkbox"/> 1	Shark drive only, left	Y00227 0
	<input type="checkbox"/> 2	Shark lights, right	Y00230 120
	<input type="checkbox"/> 2	Shark lights, left	Y00229 120
	<input type="checkbox"/> 3	Shark seat & lights, right	Y00232 175
	<input type="checkbox"/> 3	Shark seat & lights, left	Y00231 175
	<input type="checkbox"/> 4	DX G90A joystick module, right	Y00241 250
	<input type="checkbox"/> 4	DX G90A joystick module, left	Y00240 250
	<input type="checkbox"/> 5	DX 12-key flat, centre	Y00239 665
	<input type="checkbox"/> 6	Attendant control / dual joystick (not i.c.w. shark drive only), right	Y00222 440
	<input type="checkbox"/> 6	Attendant control / dual joystick (not i.c.w. shark drive only), left	Y00221 440
	<input type="checkbox"/>	<b>Fixed bracket</b>	<b>Y00264 Basic</b>
	<input type="checkbox"/> 7	Swing-away bracket	Y00265 45
	<input type="checkbox"/> 8	Mushroom knob	Y00267 42
	<input type="checkbox"/> 9	Embossed keys (not i.c.w. Shark control)	Y00263 32
	<input type="checkbox"/> 10	Push button switch small, red (only i.c.w. DX G90A)	Y00271 160
	<input type="checkbox"/> 10	Push button switch small, green (only i.c.w. DX G90A)	Y00272 160
	<input type="checkbox"/> 10	Push button switch big, red (only i.c.w. DX G90A)	Y00273 160
	<input type="checkbox"/> 10	Push button switch big, green (only i.c.w. DX G90A)	Y00274 160











**SEAT CUSHIONS**

		Option	Excl.
	<input type="checkbox"/>	SEAT CUSHION WIDTH (SW) X DEPTH (SD)	
	<input type="checkbox"/>	42 x 50 cm / 16.5 x 19.6" (SW 38 - 42 cm / 15 - 16.5", SD 44 - 52 cm / 17.3 - 20.5")	Y00784 0
	<input type="checkbox"/>	<b>46 x 50 cm / 18 x 19.7" (SW 42 - 47 cm / 16.5 - 18.5", SD 44 - 52 cm / 17.3-20.5)</b>	<b>Y00785 Basic</b>
	<input type="checkbox"/>	50 x 50 cm / 19.7 x 19.7" (SW 46 - 51 cm / 18.1 - 20.1", SD 44 - 52 cm / 17.3 - 20.5")	Y00786 0
	<input type="checkbox"/>	42 x 54 cm / 16.5 x 21.3" (SW 38 - 42 cm / 15 - 16.5", SD 48 - 56 cm / 18.9 - 22")	Y00787 95
	<input type="checkbox"/>	42 x 58 cm / 16.5 x 22.8" (SW 38 - 42 cm / 15 - 16.5", SD 52 - 60 cm / 20.5 - 23.6")	Y00788 95
	<input type="checkbox"/>	46 x 54 cm / 18.1 x 21.3" (SW 42 - 47 cm / 16.5 - 18.5", SD 48 - 56 cm / 18.9 - 22")	Y00789 95
	<input type="checkbox"/>	46 x 58 cm / 18.1 x 22.8" (SW 42 - 47 cm / 16.5 - 18.5", SD 52 - 60 cm / 20.5 - 23.6")	Y00790 95
	<input type="checkbox"/>	50 x 54 cm / 19.7 x 21.3" (SW 46 - 51 cm / 18.1 - 20.1", SD 48 - 56 cm / 18.9 - 22")	Y00791 95
	<input type="checkbox"/>	50 x 58 cm / 19.7 x 22.8" (SW 46 - 51 cm / 18.1 - 20.1", SD 52 - 60 cm / 20.5 - 23.6")	Y00792 95
	<input type="checkbox"/>	54 x 50 cm / 21.3 x 19.7" (SW 50 - 54 cm / 19.7 - 21.3", SD 44 - 52 cm / 17.3 - 20.5"), choose also frame parts for SW 54 cm / 21.3"	Y00793 95
	<input type="checkbox"/>	54 x 54 cm / 21.3 x 21.3" (SW 50 - 54 cm / 19.7 - 21.3", SD 48 - 56 cm / 18.9 - 22"), choose also frame parts for SW 54 cm / 21.3"	Y00794 95
	<input type="checkbox"/>	54 x 58 cm / 21.3 x 22.8" (SW 50 - 54 cm / 19.7 - 21.3", SD 52 - 60 cm / 20.5 - 23.6"), choose also frame parts for SW 54 cm / 21.3"	Y00795 95
	<input type="checkbox"/> 1	Frame parts for SW 54 cm / 21.3", extended armrest brackets incl.	Y00810 95
	<input type="checkbox"/> 2	<b>Standard seat with grey middle section, sellaskin upholstery</b>	<b>Y00460 Basic</b>
	<input type="checkbox"/> 3	Memory foam, stretch rubber upholstery	Y00461 140
	<input type="checkbox"/> 4	Support seat, 2 cm / 0.8" thicker, stretch rubber upholstery	Y00462 160
	<input type="checkbox"/> 5	Balance seat, 2 cm / 0.8" thicker, stretch rubber upholstery	Y00463 180
	<input type="checkbox"/> 6	Fluid, 3 cm / 1.2" thicker, stretch rubber upholstery (not i.c.w. SW54, SD54 & SD58)	Y00464 210
	<input type="checkbox"/> 7	Air seat, 3 cm / 1.2" thicker, stretch rubber upholstery (not i.c.w. SW54, SD54 & SD58)	Y00465 285
	<input type="checkbox"/> 8	Pan with depth adjustable lip, sellaskin upholstery	Y00466 30
	<input type="checkbox"/>	Adaptive or other seat, choose desired deviations with adaptations	Y00777 -50
	<input type="checkbox"/> 9	Abduction support (12 x 12 x 7 cm / 4.7 x 4.7 x 2.8") (not i.c.w. central legrest)	Y00442 110

## BACKRESTS

Option















Excl.

BACKREST CUSHION WIDTH (BW) X HEIGHT (BH)		Option	Excl.			
			<input type="checkbox"/>	42 x 52 cm / 16.5 x 20.5"	Y00796	0
			<input type="checkbox"/>	<b>46 x 52 cm / 18.1 x 20.5"</b>	<b>Y00797</b>	<b>Basic</b>
			<input type="checkbox"/>	50 x 52 cm / 19.7 x 20.5"	Y00798	0
			<input type="checkbox"/>	42 x 56 cm / 16.5 x 22"	Y00799	95
			<input type="checkbox"/>	42 x 60 cm / 16.5 x 23.6"	Y00800	95
			<input type="checkbox"/>	46 x 56 cm / 18.1 x 22"	Y00801	95
			<input type="checkbox"/>	46 x 60 cm / 18.1 x 23.6"	Y00802	95
			<input type="checkbox"/>	50 x 56 cm / 19.7 x 22"	Y00803	95
			<input type="checkbox"/>	50 x 60 cm / 19.7 x 23.6"	Y00804	95
			<input type="checkbox"/>	54 x 52 cm / 21.3 x 20.5"	Y00805	95
			<input type="checkbox"/>	54 x 56 cm / 21.3 x 22"	Y00806	95
			<input type="checkbox"/>	54 x 60 cm / 21.3 x 23.6"	Y00807	95
			<input type="checkbox"/>	<b>1 Standard backrest with grey middle section, sellaskin upholstery</b>	<b>Y00348</b>	<b>Basic</b>
			<input type="checkbox"/>	2 Backrest with lumbar support, sellaskin upholstery	Y00353	140
			<input type="checkbox"/>	3 Backrest extra soft, stretch fabric upholstery	Y00354	135
			<input type="checkbox"/>	4 Webbed backrest, stretch fabric upholstery	Y00355	160
			<input type="checkbox"/>	5 Backrest support, sellaskin upholstery	Y00356	160
			<input type="checkbox"/>	6 Backrest with adjustable straps, stretch fabric upholstery	Y00357	200
			<input type="checkbox"/>	7 Backrest with active adjustable straps, stretch fabric upholstery	Y00358	205
			<input type="checkbox"/>	Backrest with active flexible straps, backrest height 54 cm, not i.c.w. SW 54 cm, only as an adaptation	Y00778	205
			<input type="checkbox"/>	8 Thoracic support adjustable (17 x 15 x 4 cm / 6.7 x 5.9 x 1.6"), right	Y00384	70
			<input type="checkbox"/>	8 Thoracic support adjustable (17 x 15 x 4 cm / 6.7 x 5.9 x 1.6"), left	Y00383	70
			<input type="checkbox"/>	9 Shoulderrest adjustable (24 x 17 x 3,5 cm / 9.4 x 6.7 x 1.4"), right	Y00380	95
			<input type="checkbox"/>	9 Shoulderrest adjustable (24 x 17 x 3,5 cm / 9.4 x 6.7 x 1.4"), left	Y00379	95
			<input type="checkbox"/>	10 Shoulderrest not adjustable (24 x 17 x 3,5 cm / 9.4 x 6.7 x 1.4"), right	Y00382	70
			<input type="checkbox"/>	10 Shoulderrest not adjustable (24 x 17 x 3,5 cm / 9.4 x 6.7 x 1.4"), left	Y00381	70

## LEGRESTS

Option








Excl.

			<input type="checkbox"/>	1 Central legrest (i.c.w. electrical adjustment choose calfrest) (limitations in knee angle i.c.w. RWD)	Y00311	0
			<input type="checkbox"/>	2 Trapezium calfrest (36 x 25 x 5 cm / 14.2 x 10.2 x 2")	Y00294	25
			<input type="checkbox"/>	COMPLETE SET		
			<input type="checkbox"/>	<b>3 Standard swing-away legrests with calf strap</b>	<b>Y00319</b>	<b>Basic</b>
			<input type="checkbox"/>	4 Standard swing-away legrests with calf plates (23 x 18 x 4 cm / 9 x 7 x 2") (not i.c.w. SH 49 and 51 cm / 19" and 20")	Y00320	130
			<input type="checkbox"/>	5 Standard swing-away heavy duty legrests	Y00322	245
			<input type="checkbox"/>	6 Comfort swing-away legrests with calf plates (23 x 18 x 4 cm / 9 x 7 x 2")	Y00323	295
			<input type="checkbox"/>	7 Comfort swing-away heavy duty legrests	Y00324	525
			<input type="checkbox"/>	8 Electrical comfort and swing-away legrests with calf plates (23 x 18 x 4 cm / 9 x 7 x 2") (not i.c.w. Shark control)	Y00325	770
			<input type="checkbox"/>	LEGREST ACCESSORIES (not i.c.w. heavy duty legrests)		
			<input type="checkbox"/>	9 Depth adjustable calf plate bracket, right	Y00301	65
			<input type="checkbox"/>	9 Depth adjustable calf plate bracket, left	Y00300	65
			<input type="checkbox"/>	10 Kneepad (10 x 10 x 3 cm / 3.9 x 3.9 x 1.2"), right	Y00304	42
			<input type="checkbox"/>	10 Kneepad (10 x 10 x 3 cm / 3.9 x 3.9 x 1.2"), left	Y00303	42
			<input type="checkbox"/>	11 Metal heel support, right	Y00291	42
			<input type="checkbox"/>	11 Metal heel support, left	Y00290	42
			<input type="checkbox"/>	12 Leather heel support, right	Y00293	42
			<input type="checkbox"/>	12 Leather heel support, left	Y00292	42
			<input type="checkbox"/>	13 Toe strap, right	Y00338	32
			<input type="checkbox"/>	13 Toe strap, left	Y00337	32
			<input type="checkbox"/>	14 Amputation support (22 x 22 x 7 cm / 8.7 x 8.7 x 2.8"), right	Y00287	115
			<input type="checkbox"/>	14 Amputation support (22 x 22 x 7 cm / 8.7 x 8.7 x 2.8"), left	Y00286	115

## ARMRESTS

Option

Excl.

			<input type="checkbox"/>	<b>1 Armrest PU (35 x 10 x 2,5 cm / 13.8 x 3.9 x 1"), right</b>	<b>Y00067</b>	<b>Basic</b>
			<input type="checkbox"/>	<b>1 Armrest PU (35 x 10 x 2,5 cm / 13.8 x 3.9 x 1"), left</b>	<b>Y00054</b>	<b>Basic</b>
			<input type="checkbox"/>	2 Armrest padded (30 x 8 x 3 cm / 11.8 x 3.1 x 1.2"), right	Y00069	5
			<input type="checkbox"/>	2 Armrest padded (30 x 8 x 3 cm / 11.8 x 3.1 x 1.2"), left	Y00056	5
			<input type="checkbox"/>	3 Armrest padded (40 x 8 x 3 cm / 15.7 x 3.1 x 1.2"), right	Y00068	5
			<input type="checkbox"/>	3 Armrest padded (40 x 8 x 3 cm / 15.7 x 3.1 x 1.2"), left	Y00055	5
			<input type="checkbox"/>	4 Armrest P-shape padded (35 x 9 x 3 cm / 13.8 x 3.5 x 1.2"), right	Y00072	5
			<input type="checkbox"/>	4 Armrest P-shape padded (35 x 9 x 3 cm / 13.8 x 3.5 x 1.2"), left	Y00059	5
			<input type="checkbox"/>	5 Hemiplegia armrest (not i.c.w. worktop and not i.c.w. remote side), right	Y00078	170
			<input type="checkbox"/>	5 Hemiplegia armrest (not i.c.w. worktop and not i.c.w. remote side), left	Y00065	170
			<input type="checkbox"/>	6 Flip-up armrests, set (not i.c.w. thoracic support)	Y00053	225
			<input type="checkbox"/>	7 Armrests angle adjustable, set	Y00052	155

