

Order form

PUMA 6

RETAIL 2012



handicare

T: +44 (0)8450 745 945

F: +44 (0)8450 745 946

mobilityorders@handicare.co.uk

www.handicare.co.uk

Valid from	1-1-2012	Version	1.0
Dealer	Date		
Reference No.			
Delivery address	Quantity		
	PRICE LIST ID	0088	

Standard configuration exclusive VAT £ 4.995,00

- 49 cm / 19.3" seat height
- 46 cm / 18.1" seat width
- 50 cm / 19.7" seat depth
- Rear wheel drive
- Grey shrouds
- Indoor castors
- G90 controller, fixed bracket
- Compact lights and indicators
- Standard seat cushion with grey middle section
- Standard backrest cushion with grey middle section
- Standard swing-away legrests with calf strap
- Polyurethane armrests
- Polyurethane side cushions
- Approx. 60 Ah batteries
- 8 Ah charger
- Transport brackets (for hooks & carbiners)
- Seat belt
- Full range of adjustment options available: seat width - seat depth - seat angle - seat height - centre of gravity - backrest angle.

Colour options can be found in section adaptations or ask your sales representative for other available colours and prices. The standard configuration below is indicated with 'basic'. If an alternative option is required please mark clearly by ticking the associated box.

THE Y-NUMBERS WHICH ARE MENTIONED BEHIND THE OPTIONS, ARE NOT ARTICLE NUMBERS. ALWAYS USE THE Y-NUMBERS IN COMBINATION WITH THE UNIT-ID NUMBER MENTIONED ABOVE.

BASE FRAME

		Option	Excl.
SEAT HEIGHT, INCL. STANDARD CUSHION (7 CM / 2.8"), USER WEIGHT EXCL.			
<input type="checkbox"/>	49 cm / 19.3"	Y00186	Basic
<input type="checkbox"/>	51 cm / 20.1"	Y00188	0
<input type="checkbox"/>	53 cm / 20.9"	Y00192	0
<input type="checkbox"/>	55 cm / 21.7"	Y00196	0
<input type="checkbox"/>	1 Front wheel drive	Y00105	0
<input type="checkbox"/>	2 Rear wheel drive	Y00106	Basic
<input type="checkbox"/>	3 Colour shrouds grey	Y00209	Basic
<input type="checkbox"/>	4 Colour shrouds yellow	Y00210	0
<input type="checkbox"/>	Castors type outdoor	Y00095	35
<input type="checkbox"/>	Castors type indoor	Y00096	Basic
<input type="checkbox"/>	Speed 10 km/h / 6 mph	Y00215	Basic
<input type="checkbox"/>	5 Suspended castors for extra stability and comfort (only i.c.w. castors type "outdoor")	Y00219	215
<input type="checkbox"/>	Castors air	Y00097	Basic
<input type="checkbox"/>	Castors puncture proof with solid inner tube (only i.c.w. castors type "outdoor")	Y00098	80
<input type="checkbox"/>	Drive wheels air	Y00107	Basic
<input type="checkbox"/>	Drive wheels puncture proof with solid inner tube	Y00108	100
<input type="checkbox"/>	6 Telescopic Safety Device (only on FWD, obligatory i.c.w. 10 km/h / 6.2 mph)	Y00220	140



The product picture shown may deviate from the standard configuration.

ISO 7176-19 EN12184 (1999)







CLASS 3








P.O.R. = Price on request, D.O.R. = Delivery time on request

Prices in UK Pounds excluding VAT. This price list will replace all previous price list versions. We reserve the right to change (price) modifications and incorrections.










ADJUSTMENTS & LIGHTING

		Option	Excl.
	<input type="checkbox"/> 1	Electrical backrest adjustment 32°	Y00385 699
	<input type="checkbox"/> 2	Electrical tilt adjustment 40° SH + 7 cm	Y00398 650
	<input type="checkbox"/> 3	Electrical high/low adjustment 31 cm / 12.2" SH + 7 cm	Y00391 950
	<input type="checkbox"/> 4	Lighting and indicators compact version	Y00395 Basic
			

CONTROLLERS

		Option	Excl.
		DX G90A: lighting and five electrical adjustments possible	
		DX 12-key flat: lighting and five electrical adjustments possible	
	<input type="checkbox"/> 1	DX G90A joystick module, right	Y00241 Basic
	<input type="checkbox"/> 1	DX G90A joystick module, left	Y00240 0
	<input type="checkbox"/> 2	DX 12-key flat, centre	Y00239 435
	<input type="checkbox"/> 3	Attendant control / dual joystick, right	Y00222 440
	<input type="checkbox"/> 3	Attendant control / dual joystick, left	Y00221 440
	<input type="checkbox"/>	Fixed bracket	Y00264 Basic
	<input type="checkbox"/> 4	Swing-away bracket	Y00265 45
	<input type="checkbox"/> 5	Mushroom knob	Y00267 42
	<input type="checkbox"/> 6	Embossed keys	Y00263 32
	<input type="checkbox"/> 7	Push button switch small, red (only i.c.w. DX G90A)	Y00271 160
	<input type="checkbox"/> 7	Push button switch small, green (only i.c.w. DX G90A)	Y00272 160
	<input type="checkbox"/> 7	Push button switch big, red (only i.c.w. DX G90A)	Y00273 160
	<input type="checkbox"/> 7	Push button switch big, green (only i.c.w. DX G90A)	Y00274 160

SEAT CUSHIONS

		Option	Excl.
	<input type="checkbox"/>	SEAT CUSHION WIDTH (SW) X DEPTH (SD)	
	<input type="checkbox"/>	42 x 50 cm / 16.5 x 19.6" (SW 38 - 42 cm / 15 - 16.5", SD 44 - 52 cm / 17.3 - 20.5")	Y00784 0
	<input type="checkbox"/>	46 x 50 cm / 18 x 19.7" (SW 42 - 47 cm / 16.5 - 18.5", SD 44 - 52 cm / 17.3-20.5)	Y00785 Basic
	<input type="checkbox"/>	50 x 50 cm / 19.7 x 19.7" (SW 46 - 51 cm / 18.1 - 20.1", SD 44 - 52 cm / 17.3 - 20.5")	Y00786 0
	<input type="checkbox"/>	42 x 54 cm / 16.5 x 21.3" (SW 38 - 42 cm / 15 - 16.5", SD 48 - 56 cm / 18.9 - 22")	Y00787 95
	<input type="checkbox"/>	42 x 58 cm / 16.5 x 22.8" (SW 38 - 42 cm / 15 - 16.5", SD 52 - 60 cm / 20.5 - 23.6")	Y00788 95
	<input type="checkbox"/>	46 x 54 cm / 18.1 x 21.3" (SW 42 - 47 cm / 16.5 - 18.5", SD 48 - 56 cm / 18.9 - 22")	Y00789 95
	<input type="checkbox"/>	46 x 58 cm / 18.1 x 22.8" (SW 42 - 47 cm / 16.5 - 18.5", SD 52 - 60 cm / 20.5 - 23.6")	Y00790 95
	<input type="checkbox"/>	50 x 54 cm / 19.7 x 21.3" (SW 46 - 51 cm / 18.1 - 20.1", SD 48 - 56 cm / 18.9 - 22")	Y00791 95
	<input type="checkbox"/>	50 x 58 cm / 19.7 x 22.8" (SW 46 - 51 cm / 18.1 - 20.1", SD 52 - 60 cm / 20.5 - 23.6")	Y00792 95
	<input type="checkbox"/>	54 x 50 cm / 21.3 x 19.7" (SW 50 - 54 cm / 19.7 - 21.3", SD 44 - 52 cm / 17.3 - 20.5"), choose also frame parts for SW 54 cm / 21.3"	Y00793 95
	<input type="checkbox"/>	54 x 54 cm / 21.3 x 21.3" (SW 50 - 54 cm / 19.7 - 21.3", SD 48 - 56 cm / 18.9 - 22"), choose also frame parts for SW 54 cm / 21.3"	Y00794 95
	<input type="checkbox"/>	54 x 58 cm / 21.3 x 22.8" (SW 50 - 54 cm / 19.7 - 21.3", SD 52 - 60 cm / 20.5 - 23.6"), choose also frame parts for SW 54 cm / 21.3"	Y00795 95
	<input type="checkbox"/> 1	Frame parts for SW 54 cm / 21.3", extended armrest brackets incl.	Y00810 95
	<input type="checkbox"/> 2	Standard seat with grey middle section, sellaskin upholstery	Y00460 Basic
	<input type="checkbox"/> 3	Memory foam, stretch rubber upholstery	Y00461 140
	<input type="checkbox"/> 4	Support seat, 2 cm / 0.8" thicker, stretch rubber upholstery	Y00462 160
	<input type="checkbox"/> 5	Balance seat, 2 cm / 0.8" thicker, stretch rubber upholstery	Y00463 180
	<input type="checkbox"/> 6	Fluid, 3 cm / 1.2" thicker, stretch rubber upholstery (not i.c.w. SW54, SD54 & SD58)	Y00464 210
	<input type="checkbox"/> 7	Air seat, 3 cm / 1.2" thicker, stretch rubber upholstery (not i.c.w. SW54, SD54 & SD58)	Y00465 285
	<input type="checkbox"/> 8	Pan with depth adjustable lip, sellaskin upholstery	Y00466 30
	<input type="checkbox"/> 9	Abduction support (12 x 12 x 7 cm / 4.7 x 4.7 x 2.8") (not i.c.w. central legrest)	Y00442 110

BACKRESTS

		Option	Excl.
		BACKREST CUSHION WIDTH (BW) X HEIGHT (BH)	
	<input type="checkbox"/>	42 x 52 cm / 16.5 x 20.5"	Y00796 0
	<input type="checkbox"/>	46 x 52 cm / 18.1 x 20.5"	Y00797 Basic
	<input type="checkbox"/>	50 x 52 cm / 19.7 x 20.5"	Y00798 0
	<input type="checkbox"/>	42 x 56 cm / 16.5 x 22"	Y00799 95
	<input type="checkbox"/>	42 x 60 cm / 16.5 x 23.6"	Y00800 95
	<input type="checkbox"/>	46 x 56 cm / 18.1 x 22"	Y00801 95
	<input type="checkbox"/>	46 x 60 cm / 18.1 x 23.6"	Y00802 95
	<input type="checkbox"/>	50 x 56 cm / 19.7 x 22"	Y00803 95
	<input type="checkbox"/>	50 x 60 cm / 19.7 x 23.6"	Y00804 95
	<input type="checkbox"/>	54 x 52 cm / 21.3 x 20.5"	Y00805 95
	<input type="checkbox"/>	54 x 56 cm / 21.3 x 22"	Y00806 95
	<input type="checkbox"/>	54 x 60 cm / 21.3 x 23.6"	Y00807 95

